

THE FEATURES

SLOPED BASE

When only one hand is available, the sloping inner base and inconspicuous overhang help food get onto the fork or spoon more easily.



THE STOP FUNCTION

Even with one hand, the three raised edges of the chopping board ensure that nothing slips when spreading bread.



INTERNAL-CONE FUNCTION

The cone shaped inner cup makes it easy to empty the cup without having to tilt your head back.



THE STOP FUNCTION

In case of shaking, tremor or hand cramps, the inconspicuous spoon guard ensures that the food stays on the spoon.



THE SYMPTOMS



eating with one hand



shaky hands (tremor)



motor impairments of the arms, hands and fingers



visual impairment

simply clever

THE DOSE FUNCTION

In case of swallowing disorders, the inconspicuous drinking lid helps by reducing the flow of fluid. Choking is prevented and the natural swallowing process is preserved.



THE THERMAL FUNCTION

Those who need longer can still enjoy their food warm.



THE BRAKE FUNCTION

When greasing or spooning with only one hand, the non-slip ring on the bottom helps to ensure that nothing slips on the table.



THE NON-SLIP FUNCTION

In case of numbness, shaky hands or moisture, the grip cushions made of anti-slip material guarantee safe holding.




100%
MADE IN
GERMANY.
ORNAMIN.CO.UK

USE YOUR TIME
FOR FEEDING
CARE

EATING- & DRINKING AIDS -

Bestsellers

PLATE WITH SLOPED BASE




SPOON GUARD




NON-SLIP BOARDS & FOOD PREPARATION AID



MUGS WITH INTERNAL CONE



KEEP WARM PLATE & BOWL



DISCREET & THERAPEUTIC DRINKING LID



ornamin Ltd.

4200 Waterside Centre
Solihull Parkway
Birmingham Business Park
Birmingham B37 7YN

www.ornamin.co.uk

Tel +44 (0) 121 667 6279

E-Mail: team@ornamin.co.uk

Made in Germany



ORNAMIN